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 1130 1131 1230 1231

Device Description

The Strassburg Sock is a static plantar flexion prevention device. It functions as a positioning aid to assist in the treatment of plantar fasciitis. The Strassburg Sock achieves its desired effect by dorsiflexing the toes, engaging the Windlass Mechanism, and allowing the plantar fascia to heal at its maximum anatomic length.

The Strassburg Sock comes in 2 standard sizes: Regular (model#1130) and Large (model #1230). Sizing is determined by measuring around the calf at the widest point. Regular size fits up to a 16in calf circumference and a Large will fit from 16-21in around. The Strassburg Sock is available in black and white color options.

Intended Use

The Strassburg Sock is a non-sterile, reusable orthotic positioning aid intended to maintain the foot in a controlled dorsiflexed position during periods of rest or sleep. By preventing plantar flexion and maintaining a gentle, sustained stretch of the plantar fascia through activation of the windlass mechanism, the device assists in reducing biomechanical stress and supports the conservative management of plantar fasciitis and associated heel pain.

Indications

The device is indicated for the following:

- Conservative management of plantar fasciitis.
- Reduction of heel pain and morning plantar discomfort associated with plantar fascia tension.
- Maintenance of foot dorsiflexion during rest to alleviate strain on the plantar fascia.

Contraindications

- Severe peripheral vascular disease or other circulatory disorders affecting the lower limb.
- Open wounds, ulcerations, skin infections, or dermatitis in the area of contact.
- Recent foot or ankle fractures, severe deformities, or post-surgical conditions where maintaining a dorsiflexed position is contraindicated.
- Known hypersensitivity to materials used in the device.
- Severe edema or swelling of the foot or ankle that prevents proper fitting or may be aggravated by compression.

Warnings and Precautions

Warnings

- The device is intended for use during rest or sleep only. It must not be used while walking or standing, as this may result in falls or injury.
- Excessive tension on the strap may cause discomfort, circulatory restriction, or nerve compression. Adjust the tension to ensure a comfortable, gentle stretch only.
- Discontinue use immediately if pain, numbness, tingling, or skin discoloration occurs, and consult a healthcare

professional.

- The device should be used only as directed by a healthcare professional when prescribed for therapeutic use.
- Do not use a damaged or worn device, as this may compromise performance or safety.
- The device is not sterile and should not be applied to open wounds or post-surgical sites.
- Keep out of reach of children and pets.

Precautions

- Ensure the device is properly positioned and comfortably fitted before each use.
- Inspect the skin regularly for signs of irritation, redness, or pressure areas, especially during initial use.
- Gradually increase duration of wear as tolerated by the patient.
- Clean the device according to the manufacturer's instructions to maintain hygiene and functionality.
- Do not expose the device to excessive heat, moisture, or direct sunlight for prolonged periods.
- If symptoms persist or worsen despite proper use, seek medical evaluation to rule out underlying pathology.

Instructions for Use

See packaging for detailed instruction.

- Remove all packaging, unfold the sock, holding it with the O or D Ring facing away from you.
- Secure the upper strap around he calf using the attached fastner.
- Pass the strap attached to the toe of the sock through the O or D Ring.
- Relax, to allow foot to relax naturally.

Cleaning
Handwash in warm water with mild soap, rinse thoroughly and air dry.

Disposal
Dispose of device in accordance with local regulations.

Residual Risks



When used as directed, the device presents a low risk to the user. However, the following minor residual risks may still occur:

- Mild skin irritation, redness, or pressure marks due to contact with the fabric.
- Temporary discomfort or muscle soreness from the gentle stretch applied to the plantar fascia or calf muscles.
- Transient numbness or tingling if the strap is applied too tightly or worn for extended periods.
- Slight restriction of blood flow if improperly adjusted.
- Risk of slipping or falling if the user attempts to walk or stand while wearing the device.
- Rare allergic reactions to materials used in the device.

These effects are typically mild and reversible. Users should discontinue use and consult a healthcare professional if significant pain, irritation, or circulatory symptoms occur.

Reporting of Serious Incidents

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority.

Symbol	Definition
	Consult Instructions
	CE marking conformity European Union
	Identifies the authorized representative in the European Union
	Indicates device is a medical device
	Identifies catalog number of device
	Date of Manufacture with County of Origin
	Unique Device Identifier
	Does not contain latex

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